

Anova™
Your audiometer
Your way

Configurations

		Manual	Automatic (HW)	Automatic (Békésy)	Air Conduction	Bone Conduction	Masking	PC data transfer	Sound Room Microphone
Anova™ S01	Manual screening audiometer	●	-	-	●	-	-	●	○
Anova™ S02	Manual & automatic screening audiometer	●	●	○	●	-	-	●	○
Anova™ D01	Manual & automatic diagnostic audiometer	●	●	○	●	●	●	●	○

● Standard ○ Optional

Further information

- The default AC headset is the DD45. Additional headset options (DD65 v2 & IP30) are available for all configurations
- For the Anova™ D01 configuration, a standard RadioEar B71 Bone Conduction transducer is supplied
- A wireless (Bluetooth Sanibel HM-E200) printer can also be purchased separately



Amplivox Ltd, 10393 W. 70th Street, Eden Prairie, MN 55344
www.amplivox.com | 888-941-4208 | info@amplivox.com

The Amplivox policy is one of continuous development and consequently the equipment may vary in detail from the description and specification in this publication.